

Essential Theory

Objective

Determine the unknown concentration of a strong acid (**HCl**) using titration against a strong base (**NaOH**) of known concentration.

Reaction



Key Points

- **Equivalence Point:** Moles of H_3O^+ = Moles of OH^- .
- **Endpoint:** Indicated by a color change of phenolphthalein from colorless to pink.

pH Monitoring

$$\begin{aligned}\text{pH} &= -\log[\text{H}_3\text{O}^+] \\ \text{pOH} &= -\log[\text{OH}^-]\end{aligned}$$

At 25°C:

$$K_w = [\text{H}_3\text{O}^+][\text{OH}^-] = 10^{-14}, \quad \text{pH} + \text{pOH} = 14$$

Titration Curve Phases:

1. **Before Equivalence:** Acid dominates, pH increases gradually.
2. **At Equivalence:** Neutral solution, $\text{pH} = 7$.
3. **After Equivalence:** Base dominates, pH increases steeply.

Indicator

Phenolphthalein changes color at $\text{pH} \sim 9$. The indicator's endpoint closely approximates the equivalence point for strong acid/base titrations.

Procedure

Preparation

1. Rinse burette with DI water, then with NaOH.
2. Prepare 50 mL of HCl analyte in a clean Erlenmeyer flask. Add 2-3 drops of phenolphthalein.

Pilot Titration

1. Add NaOH in large increments (~ 2 mL) until a permanent pink color appears.
2. Record the approximate endpoint volume.

Full Titration

1. Add 80% of the pilot volume quickly.
2. Near the endpoint, add NaOH in 0.2 mL increments until the solution turns permanently pink.

Titration with pH Monitoring

1. Use a pH meter to record pH at intervals of added NaOH.
2. Collect data points before and after the endpoint to construct a titration curve.

Post-Experiment

- Discard waste solutions, clean all glassware, and ensure proper storage of the pH probe.

Materials

- **Glassware:** 50 mL graduated cylinder, 50 mL burette, 250 mL Erlenmeyer flasks.
- **Reagents:** 0.100 M NaOH, HCl (unknown concentration), phenolphthalein indicator.
- **Equipment:** pH meter, funnel, white paper for better visualization.

Safety Information

Hazards

- **Sodium Hydroxide (NaOH):** Corrosive, causes severe skin and eye irritation.
- **Hydrochloric Acid (HCl):** Corrosive, can cause burns and respiratory irritation.
- **Phenolphthalein:** Toxic if ingested, a strong laxative.

Precautions

- Always wear gloves, goggles, and a lab coat.
- Immediately rinse skin with cold water if contact occurs.
- Wash hands thoroughly after handling reagents.

SDS Highlights

- **NaOH:**
 - **Health Hazard:** Severe burns on skin/eyes.
 - **Handling:** Use in a well-ventilated area.
- **HCl:**
 - **Health Hazard:** Irritates respiratory tract and skin.
 - **Handling:** Avoid inhaling vapors.
- **Phenolphthalein:**
 - **Toxicity:** Ingestion can cause severe gastrointestinal distress.